



BRODA[®]
Comfort is Essential

CENTRIC
TILT SEMI-RECLINER

ESSENTIAL POSITIONING AND PEACE OF MIND

Offering basic positioning, pressure redistribution and key features for moderate levels of care, the Centric 30VT is an ideal facility chair.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE
- 3 REMOVABLE WINGS FOR UPPER LATERAL SUPPORT
- 4 HEIGHT ADJUSTABLE ARMRESTS
- 5 ADJUSTABLE SEAT TO FLOOR HEIGHT
- 6 ELEVATING LENGTH ADJUSTABLE LEGREST
- 7 REMOVABLE FLIP DOWN FOOTREST

STANDARD MODEL SHOWN

PATIENT BENEFITS



Posterior tilt improves posture, redistributes pressure and aids in fall prevention.



Tilt with semi-recline accommodates hip flexion limitations, allowing patients to be out of bed.



Optional contracture sling accommodates various levels of knee contractures.



Optional mag wheels enhance maneuverability, allowing patients to socialize and engage in their surroundings.

CAREGIVER AND FACILITY BENEFITS

- Posterior tilt enhances pressure redistribution aiding in the maintenance of skin integrity
- Tilt and recline provides ease of use for frequent re-positioning, helping improve patient outcomes
- Enhanced positioning and innovative features allow for fewer transfers, promoting safe patient handling
- Proper seating alignment decreases patient fatigue and postural deviations reducing the use of restraints
- Quality construction provides superior durability reducing equipment maintenance and replacement costs

COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue



CTS Strapping flexes for individualized fit & superior comfort.


SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Degrees of Tilt	25°
Degrees of Recline	45°
Standard Seat Width	16", 18", 20" (40cm, 45cm, 50cm)
Custom Seat Width	22", 24", 26" 28"(55cm, 60cm, 65cm, 70cm)
Standard Seat Depth	17", 18.5" (43cm, 46cm)
Custom Seat Depth	15.5", 20" (39cm, 50cm)
Back Seat Height	32" (80cm)
Back Custom Height	33.5", 35"(84cm, 88cm)
Seat Height	Adjustable 19.5", 20.5", 21.5" (49cm, 51cm, 54cm)
Front Casters	5" Front total lock (12.5cm)
Rear Casters	1X5" Total lock, 1X5" directional lock (12.5cm)
Mag Wheel Upgrade	20", 22" or 24" Front or rear mount (50cm, 55cm, 60cm)
Weight Capacity	250 lbs. (114 kg)
Standard Padding	Headrest cushion, shoulder bolsters, lower side pads, upper side pads, calf pad, sole pad, seat and back pad with cover
HCPCS	E1031

UPHOLSTERY OPTIONS



RENAISSANCE
BURGUNDY



RENAISSANCE
CHARCOAL



SILVERTEX
SAPPHIRE BLUE

FRAME OPTIONS



GREY
HAMMER

COMMON ACCESSORIES

- Mag wheel package
- Solid or clear tray
- Padded tray cover
- Thigh belt
- Contracture sling
- Padded armrests
- O2 Holder
- IV Pole
- BRODA app (additional positioning package)



HEALTHCARE DIVISION

info@brodaseating.com

p. 1-844-552-7632

Brodaseating.com

©Broda USA, INC.

8/16