

**medela**   
mother's milk,  
everyday amazing™



# Advanced Personal Double Breast Pump Starter Set

The Advanced Personal Double Breast Pump Starter Set is designed to make breast milk feeding easier. Featuring 2-Phase Expression® technology, the Advanced Personal Double Breast Pump set was designed for moms who pump several times a day and works more like a breastfeeding baby.



1. Stimulation Phase – when babies first go to breast, they suck fast and light to start milk flowing.
2. Expression Phase – after milk flow or “let-down” starts, babies breastfeed with a slower, deeper suck, bringing out more milk faster.



## Breast Pump Features

### Trusted performance

from the #1 recommended breast pump brand.\*

### Compact motor unit

in a soft case.

### Adjustable speed/vacuum control

helps mom choose her most comfortable pump setting.

### One-touch let-down button

lets mom toggle between Stimulation and Expression phases for more efficient pumping sessions.

### Simple and effective

includes a double pumping kit with 24mm PersonalFit™ Breast Shields.



## Product Specifications:

Item #	Description	Qty/ Unit	Shipper Case Dimensions (LxWxH) (inches)	Case Weight (lbs)	Qty/Case (each)
101035078	Advanced Personal Double Breast Pump Starter Set	1	9" x 8.5" x 6.75"	4.82	1

**To order or for more information, contact your Medela Sales Consultant.**

(1) Instructions for use  
 (4) Disposable Nursing Pads sample (not pictured)

**Important:** Plastic bottles and component parts become brittle when frozen and may break when dropped. Bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened, or knocked over. Take appropriate care in handling bottles and components. Do not use the breastmilk if bottles or components become damaged.



Medela LLC, 1101 Corporate Drive, McHenry, IL 60050  
 Phone: (815) 349-3980 Email: [insurance@medela.com](mailto:insurance@medela.com) [www.medela.com](http://www.medela.com)

**Content appearance may vary from pictures.**

Medela reserves the right to substitute any component or accessory with a replacement of equivalent performance.  
 \*QuintilesIMS ProVoice Survey; Nov-Oct 2017.

Medela wordmark and logo, Pump in Style, 2-Phase Expression and Invented by Medela 2-Phase Expression are registered in the U.S. Patent and Trademark Office. PersonalFit and Mother's Milk, Everyday Amazing are trademarks of Medela.

1549130 A 0318 © 2018 Medela 4551088