

Revive

SHOWER/COMMODOE
BRODA[®]

Tilt & Recline Shower Commode



Shown with optional headrest.



Anterior and posterior tilt allow exceptional maneuverability.




Dynamic positioning permits for safe and thorough bathing.



Enhances patient comfort and support with full elevating legrest.

SPECIFICATIONS

Infinitely Adjustable Tilt	up to 12° posterior, up to 10° anterior
Infinitely Adjustable Recline	up to 67°
Seat Width	18", 22"
Seat Depth	18"
Weight Capacity	400 lbs.
Frame Warranty	10 years on 16 Gauge tubular stainless steel white frame



Comfort Tension Seating[®] increases support and sitting tolerance.



Concave back accommodates and supports gluteal shelf.



Accommodates positioning needs for safe and thorough bathing.



Tilt facilitates patient transfers and ease of patient care.



Comfort Tension Seating[®] increases support and sitting tolerance.



Shown with optional headrest.

SPECIFICATIONS

Infinitely Adjustable Tilt	up to 14° posterior, up to 8° anterior
Infinitely Adjustable Recline	up to 65°
Seat Width	27", up to 33"
Seat Depth	21", 23"
Weight Capacity	800 lbs.
Frame Warranty	10 years on 16 Gauge tubular stainless steel white frame