

# Tranquille

*Auto-Locking Glider*



*Shown with tall back and neckrest cushion.*



*Gentle gliding motion calms and soothes agitation.*



*Stimulates active patients and reduces the need to wander.*



*Auto-locking feature increases safety and prevents falls.*



*Comfort Tension Seating<sup>®</sup> increases support and sitting tolerance.*

## **Patient, Facility, and Caregiver Benefits**

- Provides necessary stimulation for memory care patients.
- Increases patient safety and contentment.
- Promotes proper spinal posture and pelvic positioning.

### Contentment and Independence

#### SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	Broda Comfort Tension Seating*
Wheels	Rear transport wheel, wheels for non-patient transport*
Weight Capacity	250 lbs. (114 kg)**
Standard Padding	Seat and back pad vinyl cover
Medicare HCPCS	A9270
Lifetime Warranty	10 years on frame

\*The Tranquille Auto-Locking Glider is not a transport chair.  
\*\*Adjustments available to accommodate larger patients.

Model	Seat width	Seat depth	Seat height	Back height
Small low-back 100-10AL	18"	17"	16"	26"
Low-back 100-15AL	20"	17"	18"	26"
High-back 100-20AL	20"	17"	18"	32"

#### FRAME OPTIONS



#### COMMON ACCESSORIES

- Padded armrests
- Depth & height adjustable armrest
- Full or half tray
- Clear tray
- Contoured neck rest (high back glider only)

#### UPHOLSTERY OPTIONS

